

Product Name: c-Jun Mouse Monoclonal Antibody
Catalog #: AMM81027



Summary

Production Name	c-Jun Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IHC,ICC,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	PBS containing 0.03% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	c-Jun
Alternative Names	AP1; AP-1; c-Jun; Jun
Gene ID	3725.0
SwissProt ID	P05412.Purified recombinant fragment of human c-Jun expressed in E. Coli.

Application

Dilution Ratio	WB:1:500-1:2000,IHC:1:200-1:1000,ICC:1:200-1:1000,ELISA:1:10000
Molecular Weight	43kDa

Background

This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the

Product Name: c-Jun Mouse Monoclonal Antibody
Catalog #: AMM81027




viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.


Research Area

TGF-beta signaling pathway, MAPK signaling pathway


Image Data



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Western blot analysis using c-Jun mouse mAb against NIH/3T3 cell lysate.




Immunofluorescence analysis of HeLa cells using c-Jun mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin


Product Name: c-Jun Mouse Monoclonal Antibody
Catalog #: AMM81027



filaments have been labeled with Alexa Fluor-555 phalloidin.



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using c-Jun mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissues using c-Jun mouse mAb with DAB staining.

Note

For research use only.