
Product Name: EPN1 Mouse Monoclonal Antibody**Catalog #: AMM81420**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	60.3kDa

Antigen Information


Gene Name	EPN1
Alternative Names	Epsin 1
Gene ID	29924.0
SwissProt ID	Q9Y6I3
Immunogen	Purified recombinant fragment of human EPN1 (AA: 106-254) expressed in E. Coli.

Background


The protein encoded by this gene binds clathrin and is involved in the endocytosis of clathrin-coated vesicles. Three transcript variants encoding different isoforms have been found for this gene.

Research Area

Image Data



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of A431 cells using EPN1 mouse mAb (green) and negative control (red).