

Product Name: Rab3a Mouse Monoclonal Antibody**Catalog #: AMM81550**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	25kDa

Antigen Information


Gene Name	Rab3a
Alternative Names	Rab3a
Gene ID	5864.0
SwissProt ID	P20336
Immunogen	Purified recombinant fragment of human Rab3a (AA: 1-220) expressed in E. Coli.

Background


RAB3A (RAB3A, member RAS oncogene family) is a protein-coding gene. Diseases associated with RAB3A include choroideremia. GO annotations related to this gene include protein C-terminus binding and GTP binding. An important paralog of this gene is RAB10.

Research Area

Image Data



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HeLa cells using Rab3a mouse mAb (green) and negative control (red).