

Product Name: DHX58 Mouse Monoclonal Antibody**Catalog #: AMM81675**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	76.6kDa

Antigen Information


Gene Name	DHX58
Alternative Names	LGP2; RLR-3; D11LGP2; D11lgp2e
Gene ID	79132.0
SwissProt ID	Q5TEU4
Immunogen	Purified recombinant fragment of human DHX58 (AA: 479-678) expressed in E. Coli.

Background


DHX58 (DEXH (Asp-Glu-X-His) box polypeptide 58) is a protein-coding gene. GO annotations related to this gene include single-stranded RNA binding and helicase activity. An important paralog of this gene is IFIH1.

Research Area

Image Data



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HeLa cells using DHX58 mouse mAb (green) and negative control (red).