

Product Name: TTF1 Mouse Monoclonal Antibody**Catalog #: AMM81694**

For research use only.

Summary

| | |
|----------------------|---|
| Description | Mouse monoclonal Antibody |
| Host | Mouse |
| Application | ELISA,FC |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Modification | Unmodified |
| Isotype | Mouse IgG1 |
| Clonality | Monoclonal |
| Form | Liquid |
| Concentration | 1mg/ml |
| Storage | Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles. |
| Shipping | Ice bags |
| Buffer | Purified antibody in PBS with 0.05% sodium azide |
| Purification | Affinity Purification |

Application

| | |
|-------------------------|-------------------------------------|
| Dilution Ratio | ELISA 1:5000-1:20000,FC 1:200-1:400 |
| Molecular Weight | 103kDa |

Antigen Information

| | |
|--------------------------|---|
| Gene Name | TTF1 |
| Alternative Names | TTF-1; TTF-I |
| Gene ID | 7270.0 |
| SwissProt ID | Q15361 |
| Immunogen | Purified recombinant fragment of human TTF1 (AA: 1-150) expressed in E. Coli. |


Background

This gene encodes a transcription termination factor that is localized to the nucleolus and plays a critical role in ribosomal gene transcription. The encoded protein mediates the termination of RNA polymerase I transcription by binding to Sal box terminator elements downstream of pre-rRNA coding regions. Alternatively spliced transcript variants encoding multiple


isoforms have been observed for this gene. This gene shares the symbol/alias 'TFF1' with another gene, NK2 homeobox 1, also known as thyroid transcription factor 1, which plays a role in the regulation of thyroid-specific gene expression.

Research Area

Image Data



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using TTF1 mouse mAb (green) and negative control (red).