

Product Name: UL37 Mouse Monoclonal Antibody**Catalog #: AMM81746**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,ICC,ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	IHC 1:200-1:1000,ICC 1:50-1:250,ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	120.5kDa

Antigen Information


Gene Name	UL37
Alternative Names	UL37
Gene ID	2703358.0
SwissProt ID	Q96RL1
Immunogen	Purified recombinant fragment of Human herpesvirus UL37 (AA: 970-1119) expressed in E. Coli.

Background


UL37 complexed with large tegument protein; involved in virion morphogenesis

Research Area


Image Data




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of MCF-7 cells using UL37 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using UL37 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human esophageal cancer tissues using UL37 mouse mAb with DAB staining.