

Product Name: MUC5AC Mouse Monoclonal Antibody**Catalog #: AMM81761**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	585kDa

Antigen Information

Gene Name	MUC5AC
Alternative Names	TBM; leB; MUC5; mucin
Gene ID	4586.0
SwissProt ID	P98088
Immunogen	Purified recombinant fragment of human MUC5AC (AA: 5528-5627) expressed in E. Coli.


Background

MUC5AC (Mucin 5AC, Oligomeric Mucus/Gel-Forming) is a Protein Coding gene. Diseases associated with MUC5AC include stomach cancer and pancreatic ductal carcinoma. Among its related pathways are O-linked glycosylation and Transport to the Golgi and subsequent modification. GO annotations related to this gene include extracellular matrix structural constituent. An


important paralog of this gene is OTOG.

Research Area

Image Data



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of Hela cells using MUC5AC mouse mAb (green) and negative control (red).