
Product Name: Phospho-4E-BP1 (Ser65) Mouse Monoclonal Antibody**Catalog #: AMM81829**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	12.6kDa

Antigen Information

Gene Name	Phospho-4E-BP1 (Ser65)
Alternative Names	EIF4EBP1; BP-1; 4EBP1; 4E-BP1; PHAS-I
Gene ID	1978.0
SwissProt ID	Q13541
Immunogen	Synthesized peptide of human Phospho-4E-BP1 (Ser65).

Background


This gene encodes one member of a family of translation repressor proteins. The protein directly interacts with eukaryotic translation initiation factor 4E (eIF4E), which is a limiting component of the multisubunit complex that recruits 40S ribosomal subunits to the 5' end of mRNAs. Interaction of this protein with eIF4E inhibits complex assembly and represses translation. This

protein is phosphorylated in response to various signals including UV irradiation and insulin signaling, resulting in its dissociation from eIF4E and activation of mRNA translation.


Research Area

PI3K-Akt signaling pathway, mTOR signaling pathway


Image Data




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of Jurkat cells using Phospho-4E-BP1 (Ser65) mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human stomach cancer tissues using Phospho-4E-BP1 (Ser65) mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using Phospho-4E-BP1 (Ser65) mouse mAb with DAB staining.