
Product Name: BCL9L Mouse Monoclonal Antibody**Catalog #: AMM81933**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	IHC 1:200-1:1000,ELISA 1:5000-1:20000
Molecular Weight	157kDa

Antigen Information


Gene Name	BCL9L
Alternative Names	BCL9-2; DLNB11
Gene ID	283149.0
SwissProt ID	Q86UU0
Immunogen	Purified recombinant fragment of human BCL9L (AA: 606-751) expressed in E. Coli.

Background


BCL9L (B-Cell CLL/Lymphoma 9-Like) is a Protein Coding gene. Among its related pathways are Wnt Signaling Pathways: beta-Catenin-dependent Wnt Signaling and Signaling by Wnt. GO annotations related to this gene include transcription coactivator activity and beta-catenin binding. An important paralog of this gene is BCL9.

Research Area

Image Data



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using BCL9L mouse mAb with DAB staining.