
Product Name: CD336 Mouse Monoclonal Antibody**Catalog #: AMM82159**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	30.7kDa

Antigen Information


Gene Name	CD336
Alternative Names	NCR2; LY95; NKP44; NK-p44; dJ149M18.1
Gene ID	9436.0
SwissProt ID	O95944
Immunogen	Purified recombinant fragment of human CD336 (AA: extra 22-192) expressed in E. Coli.

Background


NCR2 is a Cytotoxicity-activating receptor that may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis.

Research Area

Image Data



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD336 mouse mAb (green) and negative control (red).