

---

**Product Name: CD370 Mouse Monoclonal Antibody****Catalog #: AMM82182**

For research use only.

**Summary**

<b>Description</b>	Mouse monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	ICC,ELISA,FC
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	Mouse IgG2b
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	ICC 1:100-1:200,ELISA 1:5000-1:20000,FC 1:200-1:400
<b>Molecular Weight</b>	27.3kDa

**Antigen Information**


<b>Gene Name</b>	CD370
<b>Alternative Names</b>	CLEC9A; DNGR1; DNGR-1; UNQ9341
<b>Gene ID</b>	283420.0
<b>SwissProt ID</b>	Q6UXN8
<b>Immunogen</b>	Purified recombinant fragment of human CD370 (AA: extra 57-241) expressed in E. Coli.

**Background**


CLEC9A is a group V C-type lectin-like receptor (CTLR) that functions as an activation receptor and is expressed on myeloid lineage cells.

## Research Area


## Image Data



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using CD370 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin.



Flow cytometric analysis of HL-60 cells using CD370 mouse mAb (green) and negative control (red).