

**Product Name: CKM Mouse Monoclonal Antibody****Catalog #: AMM82358**

For research use only.

**Summary**

<b>Description</b>	Mouse monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	ELISA
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	ELISA 1:5000-1:20000
<b>Molecular Weight</b>	43.1kDa

**Antigen Information**

<b>Gene Name</b>	CKM
<b>Alternative Names</b>	CKMM; M-CK; CPK-M
<b>Gene ID</b>	1158.0
<b>SwissProt ID</b>	P06732
<b>Immunogen</b>	Purified recombinant fragment of human CKM (AA: 1-381) expressed in E. Coli.

**Background**

The protein encoded by this gene is a cytoplasmic enzyme involved in energy homeostasis and is an important serum marker for myocardial infarction. The encoded protein reversibly catalyzes the transfer of phosphate between ATP and various phosphagens such as creatine phosphate. It acts as a homodimer in striated muscle as well as in other tissues, and as a

heterodimer with a similar brain isozyme in heart. The encoded protein is a member of the ATP:guanido phosphotransferase protein family.

## Research Area

## Image Data

