

Product Name: LILRA6 Mouse Monoclonal Antibody**Catalog #: AMM82419**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	52.4kDa

Antigen Information


Gene Name	LILRA6
Alternative Names	ILT5; ILT8; CD85b; ILT-8; LILRB3; LILRB6
Gene ID	79168.0
SwissProt ID	
Immunogen	Purified recombinant fragment of human LILRA6 (AA: 52-195) expressed in E. Coli.

Background


The protein encoded by this gene may act as receptor for class I MHC antigens.

Research Area

Image Data



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of HeLa cells using LILRA6 mouse mAb (green) and negative control (red).