
Product Name: CD267 Mouse Monoclonal Antibody**Catalog #: AMM82463**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ICC,ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	31.8kDa

Antigen Information

Gene Name	CD267
Alternative Names	CVID; RYZN; TACI; CD267; CVID2; IGAD2; TNFRSF14B
Gene ID	23495.0
SwissProt ID	O14836
Immunogen	Purified recombinant fragment of human CD267 (AA: extra(1-165)) expressed in E. Coli.


Background

The protein encoded by this gene is a lymphocyte-specific member of the tumor necrosis factor (TNF) receptor superfamily. It interacts with calcium-modulator and cyclophilin ligand (CAML). The protein induces activation of the transcription factors NFAT, AP1, and NF-kappa-B and plays a crucial role in humoral immunity by interacting with a TNF ligand. This gene is located


within the Smith-Magenis syndrome region on chromosome 17. [provided by RefSeq, Jul 2008]

Research Area


Image Data



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using CD267 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin.



Flow cytometric analysis of THP-1 cells using CD267 mouse mAb (green) and negative control (red).