
Product Name: CD45R Mouse Monoclonal Antibody**Catalog #: AMM82495**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,ICC,ELISA,FC
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:1000,ICC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	147.5kDa

Antigen Information

Gene Name	CD45R
Alternative Names	PTPRC; LCA; LY5; B220; L-CA; T200; CD45; GP180
Gene ID	5788.0
SwissProt ID	P08575
Immunogen	Purified recombinant fragment of human CD45R (AA: extra 378-577) expressed in E. Coli.

Background


The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitosis, and oncogenic transformation. This PTP contains an extracellular domain, a single transmembrane segment and two tandem intracytoplasmic

catalytic domains, and thus is classified as a receptor type PTP. This PTP has been shown to be an essential regulator of T- and B-cell antigen receptor signaling. It functions through either direct interaction with components of the antigen receptor complexes, or by activating various Src family kinases required for the antigen receptor signaling. This PTP also suppresses JAK kinases, and thus functions as a regulator of cytokine receptor signaling. Alternatively spliced transcripts variants of this gene, which encode distinct isoforms, have been reported.


Research Area

Jak-STAT signaling pathway


Image Data




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using CD45R mouse mAb against Jurkat (1), Rat Serum (2), and mouse Serum (3) cell lysate.




Immunofluorescence analysis of HeLa cells using CD45R mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.




Flow cytometric analysis of HeLa cells using CD45R mouse mAb (green) and negative control (red).



Flow cytometric analysis of Jurkat cells using CD45R mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human tonsil tissues using CD45R mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using CD45R mouse mAb with DAB staining.