
Product Name: GLUL Mouse Monoclonal Antibody**Catalog #: AMM82510**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,ELISA,FC
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	42kDa

Antigen Information

Gene Name	GLUL
Alternative Names	GS; GLNS; PIG43; PIG59
Gene ID	2752.0
SwissProt ID	P15104
Immunogen	Purified recombinant fragment of human GLUL (AA: 2-121) expressed in E. Coli.


Background

The protein encoded by this gene belongs to the glutamine synthetase family. It catalyzes the synthesis of glutamine from glutamate and ammonia in an ATP-dependent reaction. This protein plays a role in ammonia and glutamate detoxification, acid-base homeostasis, cell signaling, and cell proliferation. Glutamine is an abundant amino acid, and is important to the


biosynthesis of several amino acids, pyrimidines, and purines. Mutations in this gene are associated with congenital glutamine deficiency, and overexpression of this gene was observed in some primary liver cancer samples. There are six pseudogenes of this gene found on chromosomes 2, 5, 9, 11, and 12. Alternative splicing results in multiple transcript variants.

Research Area


Image Data




Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)




Western blot analysis using GLUL mouse mAb against Jurkat (1), and mouse liver (2) cell lysate.



Flow cytometric analysis of HeLa cells using GLUL mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cerebellar tissues using GLUL mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissues using GLUL mouse mAb with DAB staining.