

**Product Name: SALL4 Mouse Monoclonal Antibody****Catalog #: AMM82512**

For research use only.

**Summary**

<b>Description</b>	Mouse monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	WB,ELISA,FC
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
<b>Molecular Weight</b>	112kDa

**Antigen Information**


<b>Gene Name</b>	SALL4
<b>Alternative Names</b>	DRRS; HSAL4; ZNF797; dJ1112F19.1
<b>Gene ID</b>	57167.0
<b>SwissProt ID</b>	Q9UJQ4
<b>Immunogen</b>	Purified recombinant fragment of human SALL4 (954-1053) expressed in E. Coli.

**Background**


The protein encoded by this gene may be a zinc finger transcription factor. Defects in this gene are a cause of Duane-radial ray syndrome (DRRS).

## Research Area


## Image Data



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Western blot analysis using SALL4 mouse mAb against A431 (1) and NIH/3T3 (2) cell lysate.



Flow cytometric analysis of HeLa cells using SALL4 mouse mAb (green) and negative control (red).