
Product Name: SARS-CoV-2-ACE2 Mouse Monoclonal Antibody**Catalog #: AMM82615**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	ELISA 1:5000-1:20000
Molecular Weight	/

Antigen Information

Gene Name	SARS-CoV-2-ACE2
Alternative Names	SARS-CoV-2-ACE2
Gene ID	
SwissProt ID	
Immunogen	Purified recombinant fragment of SARS-CoV-2-ACE2 expressed in E. Coli.

Background**Research Area**

Image Data



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)