

---

**Product Name: D117L Mouse Monoclonal Antibody****Catalog #: AMM82668**

For research use only.

**Summary**

<b>Description</b>	Mouse monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	WB
<b>Reactivity</b>	Other
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:2000
<b>Molecular Weight</b>	/

**Antigen Information**


<b>Gene Name</b>	D117L
<b>Alternative Names</b>	p17
<b>Gene ID</b>	22220336.0
<b>SwissProt ID</b>	Q89424
<b>Immunogen</b>	Purified recombinant fragment of human D117L (AA: 1-117) expressed in E. Coli.

**Background**

structural protein p17 [ African swine fever virus ]

**Research Area**

## Image Data



Western blot analysis using D117L mAb against human