

Product Name: P22 Mouse Monoclonal Antibody**Catalog #: AMM82669**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Other
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:2000
Molecular Weight	/

Antigen Information

Gene Name	P22
Alternative Names	P22
Gene ID	22220394.0
SwissProt ID	P23169
Immunogen	Purified recombinant fragment of human P22 (AA: 29-177) expressed in E. Coli.

Background**Research Area**

Image Data



Western blot analysis using P22 mAb against human