
Product Name: P54 Mouse Monoclonal Antibody**Catalog #: AMM82670**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Other
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:2000
Molecular Weight	/

Antigen Information

Gene Name	P54
Alternative Names	E183L
Gene ID	AZI94751.1
SwissProt ID	B2MUM6
Immunogen	Purified recombinant fragment of human P54 (AA: (53-184)) expressed in E. Coli.

Background**Research Area**

Image Data



Western blot analysis using P54 mAb against human