

Product Name: UPK3B Mouse Monoclonal Antibody**Catalog #: AMM82841**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ICC,ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	ICC 1:50-1:200,ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	33.8kDa

Antigen Information

Gene Name	UPK3B
Alternative Names	P35; UP3B; UPIIIB
Gene ID	105375355.0
SwissProt ID	Q9BT76
Immunogen	Purified recombinant fragment of human UPK3B (AA: 30-180) expressed in E. Coli.


Background

UPK3B is a minor component of the apical plaques of mammalian urothelium that binds and dimerizes with uroplakin-1b (UPK1B; MIM 602380), one of the major conserved urothelium membrane proteins. The other major conserved integral membrane proteins of urothelial plaques are UPK1A (MIM 611557), UPK2 (MIM 611558), and UPK3A (MIM 611559) (Deng et al.,


2002 [PubMed 12446744]).

Research Area


Image Data



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunofluorescence analysis of HeLa cells using UPK3B mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin.



Flow cytometric analysis of THP-1 cells using UPK3B mouse mAb (green) and negative control (red).