

Product Name: XPC Mouse Monoclonal Antibody**Catalog #: AMM82884**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	106kDa

Antigen Information

Gene Name	XPC
Alternative Names	XP3; RAD4; XPCC; p125
Gene ID	7508.0
SwissProt ID	Q01831
Immunogen	Purified recombinant fragment of human XPC (AA: 32-133) expressed in mammalian.


Background

The protein encoded by this gene is a key component of the XPC complex, which plays an important role in the early steps of global genome nucleotide excision repair (NER). The encoded protein is important for damage sensing and DNA binding, and shows a preference for single-stranded DNA. Mutations in this gene or some other NER components can result in Xeroderma


pigmentosum, a rare autosomal recessive disorder characterized by increased sensitivity to sunlight with the development of carcinomas at an early age. Alternatively spliced transcript variants have been found for this gene.

Research Area


Image Data



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using XPC mouse mAb against SW480(1) ,A431(2) ,T47D(3) ,HT-29(4) ,A549 (5)and C2C12(6) cell lysate.



Flow cytometric analysis of HeLa cells using XPC mouse mAb (green) and negative control (red).