
Product Name: SPATA19 Mouse Monoclonal Antibody**Catalog #: AMM82906**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	IHC 1:200-1:1000,ELISA 1:5000-1:20000
Molecular Weight	19.1kDa

Antigen Information


Gene Name	SPATA19
Alternative Names	CT132; SPAS1; spergen1
Gene ID	219938.0
SwissProt ID	Q7Z5L4
Immunogen	Purified recombinant fragment of human SPATA19 (AA: 25-167) expressed in E. Coli.

Background


SPATA19 (Spermatogenesis Associated 19) is a Protein Coding gene. Diseases associated with SPATA19 include Multiple Mucosal Neuroma.

Research Area

Image Data



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Immunohistochemical analysis of paraffin-embedded human testis tissues using SPATA19 mouse mAb with DAB staining.