
Product Name: NES Mouse Monoclonal Antibody**Catalog #: AMM83021**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS containing 0.03% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	177kDa

Antigen Information


Gene Name	NES
Alternative Names	NESTIN
Gene ID	10763.0
SwissProt ID	P48681
Immunogen	Purified recombinant fragment of human NES (AA: 419-588) expressed in E. Coli.

Background


This gene encodes a member of the intermediate filament protein family and is expressed primarily in nerve cells.

Research Area

Image Data



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of HeLa cells using NES mouse mAb (green) and negative control (purple).