

Product Name: TRAFD1 Mouse Monoclonal Antibody**Catalog #: AMM83042**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	PBS containing 0.03% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	64.8kDa

Antigen Information


Gene Name	TRAFD1
Alternative Names	FLN29
Gene ID	10906.0
SwissProt ID	O14545
Immunogen	Purified recombinant fragment of human TRAFD1 (AA: 401-582) expressed in E. Coli.

Background


The innate immune system confers host defense against viral and microbial infection, and TRAFD1 is a negative feedback regulator that controls excessive immune responses

Research Area


Image Data



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Western blot analysis using TRAFD1 mouse mAb against HEK293 (1), Raji (2), and Jurkat (3) cell lysate.



Flow cytometric analysis of HeLa cells using TRAFD1 mouse mAb (green) and negative control (purple).