
Product Name: PIWIL4 Mouse Monoclonal Antibody**Catalog #: AMM83057**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,ELISA,FC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:1000,ELISA 1:5000-1:20000,FC 1:200-1:400
Molecular Weight	96.6kDa

Antigen Information


Gene Name	PIWIL4
Alternative Names	HIWI2; MIWI2
Gene ID	143689.0
SwissProt ID	Q7Z3Z4
Immunogen	Purified recombinant fragment of human PIWIL4 (AA: 304-434) expressed in E. Coli.

Background

PIWIL4 belongs to the Argonaute family of proteins, which function in development and maintenance of germline stem cells

Research Area

Image Data



Red: Control Antigen (100ng);