

Product Name: DNAL4 Mouse Monoclonal Antibody**Catalog #: AMM83071**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,ELISA
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	IHC 1:200-1:1000,ELISA 1:5000-1:20000
Molecular Weight	12kDa

Antigen Information


Gene Name	DNAL4
Alternative Names	PIG27
Gene ID	10126.0
SwissProt ID	O96015
Immunogen	Purified recombinant fragment of human DNAL4 (AA: 1-105) expressed in E. Coli.

Background

DNAL4 is a component of the dynein motor complex

Research Area

Image Data



Red: Control Antigen (100ng);