
Product Name: Human IgG3 Mouse Monoclonal Antibody**Catalog #: AMM83089**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:30000-1:40000
Molecular Weight	Calculated MW: 41 kDa; Observed MW: 41 kDa

Antigen Information

Gene Name	Human IgG3
Alternative Names	Human IgG3
Gene ID	
SwissProt ID	P01860
Immunogen	Recombinant human IgG3

Background**Research Area**

Image Data



Microtiter wells were coated with human immunoglobulins with various isotypes. Human IgG3 Mouse Monoclonal Antibody was used as the Human IgG3 Mouse Monoclonal Antibody and peroxidase conjugated goat anti-mouse IgG was used as the secondary antibody.