
Product Name: GBP6 Mouse Monoclonal Antibody**Catalog #: AMM83092**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:200,IP 1:100-1:200
Molecular Weight	72kDa

Antigen Information


Gene Name	GBP6
Alternative Names	GBP6
Gene ID	
SwissProt ID	Q6ZN66
Immunogen	Purified recombinant fragment of human GBP6

Background**Research Area**


Image Data



Western blot analysis using GBP6 mouse mAb against tonsil tissue lysate



Lane 1: Tonsil tissue lysate ; Lane 2: GBP6 immunoprecipitated from tonsil tissue lysate by AMM83092; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83092 can immunoprecipitate GBP6;



IHC-P analysis of esophagus tissue by GBP6 antibody (AMM83092)