
Product Name: DNAJB6 Mouse Monoclonal Antibody**Catalog #: AMM83099**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	Calculated MW: 36 kDa; Observed MW: 36 kDa

Antigen Information


Gene Name	DNAJB6
Alternative Names	DNAJB6
Gene ID	
SwissProt ID	O75190
Immunogen	Purified recombinant fragment of human DNAJB6

Background**Research Area**


Image Data



Western blot analysis using DNAJB6 mouse mAb against SH-SY5Y lysate



Lane 1: SH-SY5Y lysate ; Lane 2: DNAJB6 immunoprecipitated from SH-SY5Y lysate by AMM83099 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83099 can immunoprecipitate DNAJB6;



IHC-P analysis of testis tissue by DNAJB6 antibody (AMM83099)