
Product Name: STAMBP Mouse Monoclonal Antibody**Catalog #: AMM83104**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	Calculated MW: 48 kDa; Observed MW: 48 kDa

Antigen Information


Gene Name	STAMBP
Alternative Names	STAMBP
Gene ID	
SwissProt ID	O95630
Immunogen	Purified recombinant fragment of human STAMBP

Background**Research Area**


Image Data



Western blot analysis using STAMBP mouse mAb against Jurkat lysate



Lane 1: Jurkat lysate ; Lane 2: STAMBP immunoprecipitated from Jurkat lysate by AMM83104 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83104 can immunoprecipitate STAMBP;



IHC-P analysis of testis tissue by STAMBP antibody (AMM83104)