
Product Name: USH1C Mouse Monoclonal Antibody**Catalog #: AMM83107**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:500-1:1000
Molecular Weight	62.2kDa

Antigen Information


Gene Name	USH1C
Alternative Names	USH1C
Gene ID	
SwissProt ID	Q9Y6N9
Immunogen	Purified recombinant fragment of human USH1C

Background**Research Area**

Image Data



Western blot analysis using USH1C mouse mAb against Caco-2 lysate



IHC-P analysis of small intestine tissue by USH1C antibody (AMM83107)