
Product Name: PKD2 Mouse Monoclonal Antibody**Catalog #: AMM83108**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:500-1:1000
Molecular Weight	109.7kDa

Antigen Information


Gene Name	PKD2
Alternative Names	PKD2
Gene ID	
SwissProt ID	Q13563
Immunogen	Purified recombinant fragment of human PKD2

Background**Research Area**

Image Data



Lane 1: 15 μ g of HEK-293 lysate; Lane 2: 15 μ g of MCF-7 lysate; Lane 3: 15 μ g of HeLa lysate; Lane 4: 15 μ g of SH-SY5Y lysate; Lane 5: 15 μ g of mouse kidney lysate; Result: PKD2 Mouse Monoclonal Antibody can detect PKD2 by Western blotting;



IHC-P analysis of kidney tissue by PKD2 antibody (AMM83108)