
Product Name: ENO1 Mouse Monoclonal Antibody**Catalog #: AMM83109**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:200-1:500,IP 1:100-1:200
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 47 kDa

Antigen Information


Gene Name	ENO1
Alternative Names	ENO1
Gene ID	
SwissProt ID	P06733
Immunogen	Purified recombinant fragment of human ENO1

Background**Research Area**


Image Data



Lane 1: 15 µg of SH-SY5Y lysate; Lane 2: 15 µg of U-251MG lysate; Lane 3: 15 µg of MCF-7 lysate; Lane 4: 15 µg of A-431 lysate; Lane 5: 15 µg of HeLa lysate; Lane 6: 15 µg of THP-1 lysate; Lane 7: 15 µg of Jurkat lysate; Lane 8: 15 µg of Ramos lysate; Lane 9: 15 µg of human plasma; Result: ENO1 Mouse Monoclonal Antibody can detect human ENO1 by Western blotting.



Lane 1: 15 µg of MCF-7 lysate; Lane 2: ENO1 immunoprecipitated from MCF-7 lysate by AMM83109; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83109 can immunoprecipitate ENO1;



IHC-P analysis of colon tissue by ENO1 antibody (AMM83109)