
Product Name: CARD9 Mouse Monoclonal Antibody**Catalog #: AMM83115**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	Calculated MW: 62 kDa; Observed MW: 62 kDa

Antigen Information


Gene Name	CARD9
Alternative Names	CARD9
Gene ID	
SwissProt ID	Q9H257
Immunogen	Purified recombinant fragment of human CARD9

Background**Research Area**

Image Data



Western blot analysis using CARD9 mouse mAb against THP-1 lysate



Lane 1: THP-1 lysate ; Lane 2: CARD9 immunoprecipitated from THP-1 lysate by AMM83115 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83115 can immunoprecipitate CARD9;