

Product Name: PADI2 Mouse Monoclonal Antibody**Catalog #: AMM83118**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	PADI2
Alternative Names	PADI2
Gene ID	
SwissProt ID	Q9Y2J8
Immunogen	Purified recombinant fragment of human PADI2

Background**Research Area**

Image Data



Lane 1: PC-1 lysate