
Product Name: ZBTB17 Mouse Monoclonal Antibody**Catalog #: AMM83128**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Rat, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000
Molecular Weight	87.9kDa

Antigen Information

Gene Name	ZBTB17
Alternative Names	ZBTB17
Gene ID	
SwissProt ID	Q13105
Immunogen	Purified recombinant fragment of human ZBTB17

Background**Research Area**

Image Data



Lane 1: 15 μ g of U-937 lysate; Lane 2: 15 μ g of rat lung tissue lysate; Lane 3: 15 μ g of mouse lung tissue lysate; Result: ZBTB17 Mouse Monoclonal Antibody can detect ZBTB17 by Western blotting;