

Product Name: C7orf50 Mouse Monoclonal Antibody
Catalog #: AMM83133



Summary

Production Name	C7orf50 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IP,IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	C7orf50
Alternative Names	C7orf50
Gene ID	
SwissProt ID	Q9BRJ6.Purified recombinant fragment of human C7orf50

Application

Dilution Ratio	WB:1:500-1:1000,IP:1:100-1:200,IHC:1:500-1:1000
Molecular Weight	22.1kDa


Background

Product Name: C7orf50 Mouse Monoclonal Antibody
Catalog #: AMM83133




Research Area


Image Data



Western blot analysis using C7orf50 mouse mAb against MCF7 lysate



Lane 1: MCF7 lysate ; Lane 2: C7orf50 immunoprecipitated from MCF7 lysate by AMM83133 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83133 can immunoprecipitate C7orf50;



IHC-P analysis of kidney tissue by C7orf50 antibody (AMM83133)

Product Name: C7orf50 Mouse Monoclonal Antibody
Catalog #: AMM83133



Note

For research use only.