

Product Name: CARS1 Mouse Monoclonal Antibody
Catalog #: AMM83145



Summary

Production Name	CARS1 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IP,IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	CARS1
Alternative Names	CARS1
Gene ID	
SwissProt ID	P49589.Purified recombinant fragment of human CARS1

Application

Dilution Ratio	IP:1:100-1:200,IHC:1:100-1:200
Molecular Weight	/


Background

Product Name: CARS1 Mouse Monoclonal Antibody
Catalog #: AMM83145




Research Area

Image Data



Lane 1: HEK-293 lysate ; Lane 2: CARS1 immunoprecipitated from HEK-293 lysate by AMM83145 ; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: AMM83145 can immunoprecipitate CARS1;



IHC-P analysis of pancreas tissue by CARS1 antibody (AMM83145)

Note

For research use only.