

Product Name: TCL1B Mouse Monoclonal Antibody**Catalog #: AMM83146**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:200-1:500
Molecular Weight	/

Antigen Information

Gene Name	TCL1B
Alternative Names	TCL1B
Gene ID	
SwissProt ID	O95988
Immunogen	Purified recombinant fragment of human TCL1B

Background**Research Area**

Image Data



IHC-P analysis of testis tissue by TCL1B antibody (TCL1B Mouse Monoclonal Antibody)