

Product Name: ALDH3A1 Mouse Monoclonal Antibody
Catalog #: AMM83148



Summary

Production Name	ALDH3A1 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	ALDH3A1
Alternative Names	ALDH3A1
Gene ID	
SwissProt ID	P30838. Purified recombinant fragment of human ALDH3A1

Application

Dilution Ratio	IHC:1:500-1:1000
Molecular Weight	/


Background

Product Name: ALDH3A1 Mouse Monoclonal Antibody
Catalog #: AMM83148



Research Area

Image Data



IHC-P analysis of stomach tissue by anti-ALDH3A1 antibody (AMM83148)

Note

For research use only.