
Product Name: TXNDC5 Mouse Monoclonal Antibody**Catalog #: AMM83150**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	47.6kDa

Antigen Information


Gene Name	TXNDC5
Alternative Names	TXNDC5
Gene ID	
SwissProt ID	Q8NBS9
Immunogen	Purified recombinant fragment of human TXNDC5

Background**Research Area**


Image Data



Western blot analysis using TXNDC5 mouse mAb against Hep G2 lysate



Lane 1: Hep G2 lysate ; Lane 2: TXNDC5 immunoprecipitated from Hep G2 lysate by AMM83150 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83150 can immunoprecipitate TXNDC5;



IHC-P analysis of tonsil tissue by TXNDC5 antibody (AMM83150)