
Product Name: EBAG9 Mouse Monoclonal Antibody**Catalog #: AMM83157**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	EBAG9
Alternative Names	EBAG9
Gene ID	
SwissProt ID	O00559
Immunogen	Purified recombinant fragment of human EBAG9

Background**Research Area**

Image Data



Lane 1: T-47D lysate ; Lane 2: EBAG9 immunoprecipitated from T-47D lysate by EBAG9 Mouse Monoclonal Antibody ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: EBAG9 Mouse Monoclonal Antibody can immunoprecipitate EBAG9;