

Product Name: TGM2 Mouse Monoclonal Antibody**Catalog #: AMM83158**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:500-1:1000,IP 1:100-1:200
Molecular Weight	90kDa

Antigen Information


Gene Name	TGM2
Alternative Names	TGM2
Gene ID	
SwissProt ID	P21980
Immunogen	Purified recombinant fragment of human TGM2

Background**Research Area**


Image Data



Western blot analysis using TGM2 mouse mAb against A549 lysate



Lane 1: A549 lysate ; Lane 2: TGM2 immunoprecipitated from A549 lysate by AMM83158 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83158 can immunoprecipitate TGM2;



IHC-P analysis of colon tissue by TGM2 antibody (AMM83158)