
Product Name: MYO6 Mouse Monoclonal Antibody**Catalog #: AMM83162**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	149kDa

Antigen Information


Gene Name	MYO6
Alternative Names	MYO6
Gene ID	
SwissProt ID	Q9UM54
Immunogen	Purified recombinant fragment of human MYO6

Background**Research Area**


Image Data



Western blot analysis using MYO6 mouse mAb against HT-29 lysate



Lane 1: HT-29 lysate ; Lane 2: MYO6 immunoprecipitated from HT-29 lysate by AMM83162 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83162 can immunoprecipitate MYO6;



IHC-P analysis of kidney tissue by MYO6 antibody (AMM83162)