
Product Name: ABRAXAS2 Mouse Monoclonal Antibody**Catalog #: AMM83164**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	/

Antigen Information


Gene Name	ABRAXAS2
Alternative Names	ABRAXAS2
Gene ID	
SwissProt ID	Q15018
Immunogen	Purified recombinant fragment of human ABRAXAS2

Background**Research Area**

Image Data



Lane 1: Hep G2 lysate ; Lane 2: ABRAXAS2 immunoprecipitated from Hep G2 lysate by ABRAXAS2 Mouse Monoclonal Antibody ; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: ABRAXAS2 Mouse Monoclonal Antibody can immunoprecipitate ABRAXAS2;



IHC-P analysis of colon cancer tissue by ABRAXAS2 antibody (AMM83164)