
Product Name: CCDC137 Mouse Monoclonal Antibody**Catalog #: AMM83168**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgM
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:50-1:200
Molecular Weight	/

Antigen Information

Gene Name	CCDC137
Alternative Names	CCDC137
Gene ID	
SwissProt ID	Q6PK04
Immunogen	Purified recombinant fragment of human CCDC137

Background**Research Area**

Image Data



IHC-P analysis of colon cancer tissue by CCDC137 antibody (CCDC137 Mouse Monoclonal Antibody)