

---

**Product Name: SCAMP5 Mouse Monoclonal Antibody****Catalog #: AMM83174**

For research use only.

**Summary**

<b>Description</b>	Mouse monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	WB,IHC,IP
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
<b>Molecular Weight</b>	26.1kDa

**Antigen Information**


<b>Gene Name</b>	SCAMP5
<b>Alternative Names</b>	SCAMP5
<b>Gene ID</b>	
<b>SwissProt ID</b>	Q8TAC9
<b>Immunogen</b>	Purified recombinant fragment of human SCAMP5

**Background****Research Area**


## Image Data



Western blot analysis using SCAMP5 mouse mAb against mouse brain lysate



Lane 1: Mouse brain lysate ; Lane 2: SCAMP5 immunoprecipitated from mouse brain lysate by AMM83174 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83174 can immunoprecipitate SCAMP5;



IHC-P analysis of cerebellum tissue by SCAMP5 antibody (AMM83174)