
Product Name: VEGFC Mouse Monoclonal Antibody**Catalog #: AMM83177**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human, Rat, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 42 kDa

Antigen Information


Gene Name	VEGFC
Alternative Names	VEGFC
Gene ID	
SwissProt ID	P49767
Immunogen	Purified recombinant fragment of human VEGFC

Background**Research Area**

Image Data



Western blot analysis using VEGFC mouse mAb against rat heart tissue lysate



Lane 1: Rat heart lysate ; Lane 2: VEGFC immunoprecipitated from rat heart lysate by AMM83177 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83177 can immunoprecipitate VEGFC;